


**Dentium Instruments**  
for Total Solution

Catalog & Manual

# Help Kit

- Easy solution for critical problems which may occur in the prosthetic process consist of 5 tools in a kit
- (Screw Remover/ Abutment Hex Remover/ Screw Tap Repair / Fixture Remover / Cover & Abutment Screw Remover)
- Compatible with most dental implant products now available on the global market
- Heavy duty with robust design and proven materials

*SuperLine*



XIH

## Screw Remover

L	Art. No.
25	XRF S2S
35	XRF S2



## Abutment Remover

L	Art. No.
20	XRHR 20
25	XRHR 25



## Screw Tap Repair

Type	Art. No.
Tap	XRSTR
11° Guide	XRSG11
8° Guide	XRSG8




## Fixture Remover

Type	Art. No.
Remover	XRFRTF
	XRFRT
	XRFRT0
Wrench	XRFRW




## Cover & Abutment Screw Remover

L	Art. No.
25	XRRHD



# NR Line



GXIH

## Screw Remover

L	Art. No.
29	GSRS
33	GSRL




## Abutment Remover

L	Art. No.
20	GAHS
25	GAHL





## Screw Tap Repair

L	Art. No.
25	GSTRP
GUIDE	GSTG10



## Fixture Remover

Type	Art. No.
25	GFXRS
30	GFXRL
Wrench	XRFRW



## Cover & Abutment Screw Remover

L	Art. No.
25	GCASR



## Ratchet adapter

L	Art. No.
13.9	GRA



## Screw Remover

### Application

To remove the remaining screw when the abutment screw is broken inside the fixture

### Advantage


Easy to remove the broken screw, as well as protect the internal threads of the fixture from being damaged

### Usage

1. Set the torque of the implant motor to 30~50 rpm in a CCW (counterclockwise) direction
2. Assemble the tool with the hand-piece
3. Run the motor while keeping the tip of the tool appropriately contacted with the broken screw until successfully removed


\*Caution: Do not overload the tool with pressure; apply moderate pressure

1



Dual Abutment


2



Use the friction force of the tool rotating counterclockwise to remove the screw

Hand-piece torque: 30~50rpm / Reverse

3



Allow the screw to gradually come out in a swaying motion

4



## Abutment Hex Remover

### Application

To remove the remaining hex when the hex portion of an abutment is broken

### Advantage


Easy to remove the broken hex, as well as protect the internal threads of the fixture from being damaged

### Usage

1. Insert the tool inside into the remaining hex hole of the fixture inside
2. Assemble the ratchet with the tool and rotate it in a CW (clockwise) direction to lock the tool tip with the remaining hex
3. Disengage the ratchet and remove the remaining hex by gently rocking the tool
4. If necessary, the hole located in the upper portion of the tool may be used with the crown ejector (not included)


\*Caution: Do not overload the tool with pressure; apply moderate pressure

1




Dual Abutment (Hex)

2



Rotate the tool clockwise so that the remaining hex gets tightly engaged to the tool

3



Once the tool is tightly locked to the hex remnant, disengage the ratchet  
Gently rock the tool until the hex is successfully removed

4



## Screw Tap Repair

### Application

To recreate the internal thread lines of the fixture when it is damaged

### Advantage

Easy to recreate the internal threads with the help of the guides corresponding to different internal angulation (8, 11 degrees) of the fixture

### Usage

1. Place the guide with corresponding degree to the fixture
2. Assemble the tap tool with ratchet
3. Start tapping using the tap tool with appropriate torque
4. If excessive debris accumulates, pause tapping and remove using suction
5. Repeat steps 3 and 4 until completed

\*Caution: Do not apply excessive torque onto the tap tool

It is highly recommended to use the ratchet after the initial engagement of the tool and the internal threads

1



IMPLANTIUM / SuperLine 11°



SimpleLine II 8°

3



Tap with the guide attached

4



Remove the tool and the guide to suction the debris

\*If excessive debris accumulates, pause tapping and remove using suction

## Fixture Remover

### Application

To remove the fixture when critically damaged with no other recovery options

### Advantage

Easy to remove the failed fixture without causing damage to the adjacent bone

### Usage

1. Assemble the tool with ratchet, and insert it into the failed fixture to be removed
2. Gently rotate the ratchet in a CCW direction until the tool is tightly locked into the fixture
3. Continue to rotate the ratchet with greater torque in a CCW direction until the failed fixture is completely removed
4. Separate the tool from the removed fixture by rotating it in a CW direction. If necessary, use the wrench (included) to hold the fixture while rotating the tool with ratchet in a CW direction

\*Caution: Sufficient irrigation should be applied to the tool to prevent excessive heating during the procedure

1 Art No. XRFRT



IMPLANTIUM / SuperLine 11°

Art No. XRFRT0



SimpleLine II 8°

2~3



Rotate the tool in a counter clockwise direction until it is tightly locked into the fixture. Continue to rotate with additional torque until the failed fixture is completely removed

4



Separate the tool from the fixture using the ratchet and the wrench that are included in the kit

## Cover & Abutment Screw Remover

### Application

To disengage the cover screw, healing abutment and abutment screw from the fixture when the 1.28 hex on the head is stripped or damaged

### Advantage

Easy to disengage the cover screw, healing abutment and abutment screw with stripped or damaged hex


### Usage

1. Assemble the tool with the ratchet and place it over the damaged 1.28 hex of the cover screw, healing abutment or abutment screw that needs to be removed
2. Gently rotate the ratchet in a CCW direction to tightly engage the tapered top of the tool into the damaged 1.28 hex.
3. Continue to rotate the ratchet in a CCW direction with greater torque until the cover screw, healing abutment or abutment screw is completely removed
4. After the removal, rotate the ratchet in a CW to separate the tool and the removed component

1 Cover Screw




2 Loading downward




Rotate the tool counterclockwise until tightly locked into the 1.28 hex of the cover screw

3 Abutment Screw



4 Loading downward



Rotate the tool counterclockwise until tightly locked into the 1.28 hex of the abutment screw

# DENTIUM LONG-TERM CLINICAL DATA

2002

2003


2004

2005

2006

2007

2008



2002. 05. 17  
Pre-op



2002. 09. 04  
Post-op



2003. 03. 15  
Final prosthesis

# Dentium

For Dentists By Dentists

2009

2010

2011

2012

2013


2014

2015

11 YEARS



2008. 04. 14  
5 years



2013. 12. 05  
11 years

over  
**10** years  
of Long  
term  
data

OVER A **DECADE** OF  
COMMITMENT TO  
THE **BEST PRODUCTS**  
FOR DENTISTS AND  
PATIENTS

**Dentium**  
For Dentists By Dentists

**Dentium**  
For Dentists By Dentists

# Dentium Instruments

for Total Solution

Catalog & Manual

**Dentium**  
For Dentists By Dentists

Specifications are subject to change without any notice.  
Some products listed in this catalog are not available in the market due to pending approval.

---

HEAD OFFICE

6F, 29, Teheran-ro 87-gil, Gangnam-gu, Seoul, Korea (135-879) Tel +82-2-555-3750 Fax +82-2-501-9560

HOME PAGE

[www.dentium.co.kr](http://www.dentium.co.kr)

INCO-1603 [Rev.1]