




It's more than portable. It's pocketable!

K2 Mobile™




Features


- Stable output power by HECOSA™ Technology
3.5W Continuous wave 6W Pulse mode



- Intra-oral illumination
- Bright & clean OLED display
- Gravity sensor for both hands
- Ideal weight and size for handheld device -19cm, 108g
- Long lasting battery life -1700mAh Lithium-ion
- Various reusable tips (auto clvable) for endo, perio, surgery Tx, etc
- Least consumables
- Adjustable fiber length



- Head cap detachable and auto clvable



- World first movable fiber for least consumables and sumables and versatile usages



- No powercord and no footpedal
- Soft button technology
- Auto emission button



- Easy menu navigation



- Voice confirmation
- User interactive



- Fatigue-free and Ergonomic design
(2014/2015 Good Design Award Winning)



- Fiber cleaver inside Cleaver blade to cut optic fiber for least consumables

