



Silfradent Surgybone Piezoelectric
Ultrasonic Bone Surgery System

SURGYBONE




silfradent



SURGICAL HANDPIECE

- Micrometric cutting
- Selective cutting
- Cavitation effect
- Minimum tissue overheating
- High visibility in the operating room
- Protection of soft tissues
- Reduction of pain
- Early rehabilitation
- High yield ceramic disks
- External irrigation set
- Light and noiseless
- Manufactured in tested material for more than 20.000 sterilization cycles

TIP OPTION



- 1 Curette
- 2 Saw
- 3 PERIONTAL SURGERY
- 4 HPISE 1
- 5 HPISE 2
- 6 HPISE 3

More Safer! More Faster! Hydric sinus lift with Water pressure

Special Tip from Pf. Sohn



- 1) Bone scrapig & Bone grafting
- 2) Ridge split & Wisdom tooorh Extraction
- 3) Periodontal surgery
- 4) Sinus lifting 1st
- 5) Sinus lifting 2st & Lateral window opening (Stopper)
- 6) Crestal approach sinus lifting



High power handpiece with internal spray



Pump with easy insertion of the irrigation set



Dynamometric wrench



Surgical tray



Kit with 6 insert tips (5 sharp, 1 smoothing)

TECNICAL DATA



- Power source: 230 V - 50/60 Hz:
- Nominal power consumption: 170 VA
- Maximum vibration: 200 micron
- Hand piece cable length: 2000 mm
- Ultrasound frequency: 22 - 38 KHz
- Footswitch: pneumatico
- Hydraulic circuit capacity: 0 - 50 ml/min
- Weight: 4,9 Kg
- Dimensions: L x P x H 270x235x170 mm

1 year

WARRANTY
FROM SHIPPING DATE



RETURN
GUARANTEE

www.minimaximplant.com.au

www.minimaxeshop.com