




VSP



EXPAND YOUR DIGITAL WORKFLOW



LIGHT, COMPACT AND USER-FRIENDLY SCANNER



PROVIDE EASY INTERFACE WITH A 4-INCH TOUCH SCREEN



IMAGE PREVIEW ON A 4-INCH TOUCH SCREEN



SLEEP MODE SUPPORT FOR POWER SAVING



SUPPORT OF MULTI-USER PRACTICES



SUPPORT STAND-ALONE MODE WITH BUILT-IN MEMORY



IMAGES CAN BE SENT TO A DESIGNATED PC ON THE NETWORK LIST



IMAGES ARE SAVED ON A BUILT-IN MEMORY CARD UP TO 500 IMAGES



EASY TO INSTALL AND USE (NETWORK SET UP AUTOMATICALLY)



HIGH DURABILITY MINIMIZES AFTER-SALES SERVICE OCCURRENCES

A SMALL, LIGHT AND COMPACTLY DESIGNED SCANNER

Vatech VSP is small and light so it can be set up in a small-sized dental clinic efficiently.

As a high-end device, Vatech VSP, has a 4-inch LCD touch screen which provides the largest image preview. The touchscreen offers quick and easy access to features.

Moreover, the ethernet cable enables rapid image transfer from scanner to pc. We developed Vatech VSP as a simple, fast, and easy way to diagnose a patient on chairside.



SPECIFICATIONS (Vatech VSP)

Supported IP	Size 0 / 1 / 2 / 3
Resolution (lp/mm)	Super Resolution: 20lp/mm, High Resolution 10lp/mm
Display	4-inch touch screen
Interface	Ethernet
Multi-user practices	Yes
Sleep mode	Yes
Weight	3.3kg
Dimensions (HxWxDmm)	233 x 141 x 191.5