



# OK-STRIPS

## Ortho K Strips



ISO13485 CE2195 FDA

# OK-STRIPS

OK strips are attached to the EVA type handpiece to make interdental removal as much as desired on the desired area.



interproximal stripping for Clear Aligner or normal dental treatments.

Non Hole Type			Hole Type		
Orthok-Strips 15SW	Orthok-Strips 15	Orthok-Strips 25	Orthok-Strips 40	Orthok-Strips 60	Orthok-Strips 80
15sw/ $\mu$ m	15/ $\mu$ m	25/ $\mu$ m	40/ $\mu$ m	60/ $\mu$ m	80/ $\mu$ m
Thickness	Thickness	Thickness	Thickness	Thickness	Thickness
0.09mm	0.09mm	0.11mm	0.13mm	0.21mm	0.27mm

## How to use



※ **Precautions** - Since dental handpieces and OK strips are very sensitive to external shocks, it is recommended not to apply excessive force and use the handpiece speed at less than 15,000rpm (medium).

## OK-STRIPS KIT

### HandPiece Set



- HandPiece
- Pin
- Strip Handle
- Space Gauge
- Strip case(4hole)
- OK-Strips(2pcs)
- OK-Strips(4pcs)\*2

## OK-STRIPS ITEM

### 1.General Type(Non hole)

Double side	Single side	Coating particle thickness
OK-S-15	OK-S-15S	15μm
OK-S-25	OK-S-25S	25μm
OK-S-40	OK-S-40S	40μm
OK-S-60	OK-S-60S	60μm
OK-S-80	OK-S-80S	80μm

### 2.Hole Type

Double side	Coating particle thickness
OK-S-15H	15μm
OK-S-25H	25μm
OK-S-40H	40μm
OK-S-60H	60μm
OK-S-80H	80μm

### 3.Saw Blade Type

Double side	Single side	Coating particle thickness
OK-S-15SW	N/A	15μm

# OK-PROXFILE



**Prosthodontic case**

Resin overhanging removal

**Orthodontic case**

Remove residual resin after bracket removal, IPR, tooth trimming & shaping

## PROXOFIELD TYPE

Double side	Single side	Coating particle thickness
OK-S-P15	N/A	15µm
OK-S-P40	N/A	40µm

**How to use** (Remove residual resin after bracket removal case)



15605 경기도 안산시 단원구 별망로 419(목내동) 광양프론티어밸리 1차 302호

302,419, Byeolmang-ro, Danwon-gu, Ansan-si, Gyeonggi-do15605, Republic of Korea

**TEL.** 82+2.6381.8160~1 **FAX.** 82+2.772.9095 **E-MAIL.** jaintek9091@gmail.com / jain98@naver.com