


# ShapeShifter™

DENTAL BLADES THAT CHANGE YOUR OUTCOMES




Surgistar

# ShapeShifter™

DENTAL BLADES THAT CHANGE YOUR OUTCOMES




Specially designed for use in Endodontic Microsurgery.



Malleable shank can be contoured to your preference. Negotiate turns and difficult angles with ease.



Narrow blades are 60% smaller than Standard, improving visibility under a microscope and allowing for more precise incisions.









Size comparison of Narrow, Standard and traditional scalpel.


# Surgistar

Wide variety of profiles and types available,  
including Standard and Narrow, both  
Malleable and Non-Malleable.


## Malleable Standard Profile

|  | Screw-In Blade | Description                            | Pkg Contents |
|--|----------------|--|--------------|
|  | S6300          | Spear, Straight                        | 6/Box        |
|  | S6400          | Round, Sharp One Side, Straight        | 12/Box       |
|  | S6500          | Lance 15°, Straight                    | 12/Box       |
|  | S6700          | Pointed, Straight                      | 12/Box       |
|  | S6900          | Round, Sharp All Around, Straight      | 6/Box        |
|  | S6915          | Round, Sharp All Around,<br>Angled 10° | 6/Box        |

## Malleable Narrow Profile

|  | Screw-In Blade | Description                             | Pkg Contents |
|--|----------------|---|--------------|
|    | USM6300        | Spear, Straight                         | 6/Box        |
|  | USM6310        | Spear, Angled 10°                       | 6/Box        |
|    | USM6400        | Round, Sharp One Side,<br>Straight      | 6/Box        |
|  | USM6700        | Pointed, Straight                       | 6/Box        |
|  | N6900          | Nordland, Sharp All Around,<br>Straight | 6/Box        |
|  | USM6900        | Round, Sharp All Around, Straight       | 6/Box        |
|  | USM6910        | Round, Sharp All Around,<br>Angled 10°  | 6/Box        |

## Malleable Mirrors

|  | Screw-In Blade | Description            | Pkg Contents |
|--|----------------|------------------------|--------------|
|  | MR2            | Round, 3.8mm           | 1/Bag        |
|  | MR1            | Round, 5mm             | 1/Bag        |
|  | MR4            | Straight Sides, Small  | 1/Bag        |
|  | MR3            | Straight Sides, Medium | 1/Bag        |
|  | MR6            | Straight Sides, Large  | 1/Bag        |






## Handles

|   | Screw-In Blade | Description                   | Pkg Contents |
|---|----------------|-------------------------------|--------------|
|  | S1301          | Blade Wrench                  | 1/Bag        |
|   | S1305          | Collet, Universal Cross-Cut   | 1/Bag        |
|    | S1310          | Round, 10cm                   | 1/Bag        |
|   | S1350          | Round, 13cm                   | 1/Bag        |
|   | S1320          | Round, 15cm                   | 1/Bag        |
|   | S1360          | Round, 5" Long, 1/4" Diameter | 1/Bag        |
|    | S1330          | Hex, 10cm                     | 1/Bag        |
|   | S1340          | Hex, 15cm                     | 1/Bag        |

## Non-Malleable Narrow Profile

|  | Screw-In Blade | Description                         | Pkg Contents |
|--|----------------|-------------------------------------|--------------|
|  | 6961           | Round, Sharp All Around, Straight   | 6/Box        |
|  | 6962           | Round, Sharp All Around, Angled 25° | 6/Box        |
|  | 6966           | Round, Sharp All Around, Angled 55° | 6/Box        |

## Non-Malleable Standard Profile

|  | Screw-In Blade | Description                         | Pkg Contents |
|--|----------------|-------------------------------------|--------------|
|    | 6300           | Spear, Straight                     | 6/Box        |
|    | 6400           | Round, Sharp One Side, Straight     | 12/Box       |
|  | 6500           | Lance 15°, Straight                 | 12/Box       |
|  | 6700           | Pointed, Straight                   | 12/Box       |
|  | 6900           | Round, Sharp All Around, Straight   | 6/Box        |
|  | 6910           | Round, Sharp All Around, Angled 10° | 6/Box        |

Distributed By:

To place an order, check on delivery status, request a sample or other customer service inquiries:  
 2310 La Mirada Dr., Vista, CA 92081  
 Tel: 1-760-598-2480 Fax: 1-760-598-2481  
 customerservice@surgistar.com



Ask your distributor for other Surgistar products!

# Surgistar

*The Cutting Edge of Tomorrow*

[www.surgistar.com](http://www.surgistar.com)